Developed By P-I Brånemark

## INSTRUCTIONS FOR USE ENGLISH

P-I Dental Components

This Instructions for Use applies to the below listed P-I Dental Components, sterile and non-sterile, included in the P-I Dental Implant Systems.

Interface	Materials	⊖ HEX	⊙ AMP	⦿ MT	Single Use	Sterile
<b>Cover Screws</b>	Titanium	G4	G4	G4	◦	◦
<b>Healing Abutments</b>	Titanium	G4	G4	G4	◦	◦
<b>Abutments</b>						
Conical Straight	Titanium	G4	G5	G5	◦	◦
Conical Angled	Titanium	G5	G5	G5	◦	◦
Cemented Cylinder	Titanium	G4	G5	G5	◦	◦
Esthetic	Titanium	G4	G5	G5	◦	◦
Contour	Titanium			G5	◦	◦
<b>Cylinders over Implant</b>	Titanium	G4	G5	G5	◦	◦
	CoCrMo	◦	◦	◦	◦	◦
Castable	POM	◦			◦	
<b>Links</b>	Titanium		G5	G5	◦	◦
<b>Scan Bodies</b>	PEEK	◦	◦	◦	◦	
<b>Healing Caps</b>	POM	◦	◦	◦	◦	
<b>Cylinders over Abutment</b>	Titanium	G4	G4	G4	◦	◦
	CoCrMo	◦	◦	◦	◦	◦
Castable	POM	◦	◦	◦	◦	
<b>Impression Copings</b>	Titanium	G5	G5	G5	◦	◦
	POM				◦	◦
<b>Screws</b>	Titanium	G5	G5	G5	◦	◦
<b>Analogs</b>	Titanium	G5	G5	G5	◦	

### Materials

G4 | commercially pure titanium grade 4 (ASTM F67).

G5 | titanium grade 5 (ASTM F136).

CoCrMo | Cobalt Chromium Molybdenum (ASTMF1537-11).

POM | Polyacetal.

PEEK | Polyether ether ketone.

Some devices may feature carbon coating.

The P-I Dental Components have various dimensional and geometric options. Anodized prosthetic Components do not indicate dimension or compatibility with other Components or Implants.

Consult the device label and options at [www.pibranemark.com](http://www.pibranemark.com)

### DISCLAIMER OF LIABILITY

Some conditions, whether combined or not, represent contraindications, limitations and risks (relative and absolute) for the treatment of patients with implants. The procedures for placement of implants are complex and require specialized training which lies in the responsibility of the respective qualified professional. The professional is the sole responsible to determine if the treatment, procedures and devices are suitable for the specific clinical case and patient. P-I disclaims all and any liability and responsibility with regards to the professional decisions and procedures when using the P-I devices.

Some devices may not be available or cleared for clinical use in your region.

### INDICATIONS FOR USE

The P-I Dental Components are part of the P-I Dental Implant Systems and are intended to be used for anchorage of crowns, bridges or dentures to the surrounding bone in the upper and lower jaws to reestablish esthetic and masticatory function. The procedures may range from prosthetic replacement of a single tooth to entire arch of bridgework. Crowns and bridges can be screw-retained and/or cemented to implants and/or abutments.

### INDICATIONS

**Cover Screws** are temporarily used on the platforms of the respective Implant during the submerged healing phase to prevent contamination, entry of residues and tissue inside the Implants.

**Healing Abutments** are temporarily used on the platforms of the respective Implant, with diameter and height options, parallel or divergent (concave) emergence, orienting and conditioning soft tissue while providing access to the respective Implant platform when using a 1 or 2-stage surgical technique.

Straight and Angled **Conical Abutments** are used on the platforms of the respective Implant interfaces, with height options to support multiple screw-retained prostheses. The maximum occlusal angulation between two Conical Abutments is 40°. The mesiodistal and interocclusal space requirements, from the Conical Abutment platform, adding the Cylinder Conical Abutment, infrastructure and esthetic material should be considered. Straight Conical Abutments associated with engaging Components (ENG) can be used for single unit screw-retained prosthesis on the AMP and MT interfaces.

**Abutment Cemented Cylinder, Esthetic Abutment** and **Contour Abutment** are used on the platforms of the respective Implant interfaces, with diameter and height options, to support cement-retained prosthesis. Esthetic and Contour Abutments feature contoured emergence and angulation options. Contour Abutments are only available for the MT interface. The Abutment Cemented Cylinder "0 height", without margin, available for the AMP and MT interfaces, is not compatible with the healing caps, impression copings, analogs and respective castable cylinders.

**Cylinders over Implant in Titanium** (Ti) are used on the platforms of the respective Implant interfaces, with diameter options, engaging (ENG), to support prostheses. Cylinders over Implant in titanium can be adjusted and feature structures for retention of esthetic material.

**Cylinders over Implant in Cobalt Chromium Molybdenum** (CoCrMo) have an upper castable part in Polyacetal – POM and are used on the platforms of the respective Implant interfaces, with diameter options, engaging (ENG), to support prostheses. The Cylinders over Implant feature complete structures for the retention of esthetic material, and the superstructure should be constructed by the laboratory.

**Links** are used on the platforms of the respective Implant interfaces, with height options, to support prostheses designed and manufactured by a laboratory using CAD/CAM systems.

**Scan Bodies** for Implants or Abutments are used on the platforms of the respective Implant or Abutment, to transfer the correct position of the Implants or Abutments using scanners to support prosthesis manufactured by a laboratory using CAD/CAM.

**Healing Caps** are used on the respective Abutments, orienting the healing and conditioning of soft tissue, preceding the installation of subsequent Components.

**Cylinders - Conical Abutment in Titanium** (Ti) are used on the respective Conical Abutments platforms to support prostheses. Cylinders Conical Abutment in titanium can be prepared and feature complete structures for esthetic material retention. The option of engaging Cylinders Conical Abutment (ENG) is only available for straight Conical Abutments with the AMP and MT interfaces.

**Cylinders - Conical Abutment in Cobalt Chromium Molybdenum** (CoCrMo) have an upper castable part in Polyacetal – POM and are used on the platforms of the respective Conical Abutments to support prostheses. The Cylinders Conical Abutment feature complete structures for the retention of esthetic material and the infrastructure should be constructed by the laboratory. The option of engaging Cylinders Conical Abutment (ENG) is only available for straight Conical Abutments with the AMP and MT interfaces.

**Impression Copings**, Open Tray (OT) or Closed Tray (CT), are used on the platforms of the respective Implants or Abutments platforms, to transfer the correct position of Implants or Abutments using impression procedures to construct a laboratory model.

**Screws** supplied separately are used for replacement, or in laboratorial procedures, to fixate the respective Component to the corresponding Implant or Component. All P-I Components are supplied with the applicable Screws.

**Analogs** and **Analogs + Digital** replicate the respective Implant or Component platform to support laboratorial procedures including the construction of models.

**Castable** Components in Polyacetal (POM) are used on the respective platform of Components or Implant after modification by a laboratory.

### CONTRAINDICATIONS

Pre-operative patient evaluation is necessary to determine any factors which may put the patient at risk from the implant placement procedure itself, or factors that may affect healing capacities of either the bone or associated soft tissue.

Components and Implants should not be used:

- In patients with a health history that contraindicates surgical intervention
- Where bone is of insufficient quality and/or quantity or not available to produce adequate stability and support
- When there is risk for overload, unfavorable jaw relations and/or parafunction
- Where there is no space to place sufficient number of implants in an optimum position to support biomechanical loads
- Where bone quantity and/or quality is not able to provide adequate implant initial stability
- Where there is a proven allergy or hypersensitivity to titanium, one of the metals in the titanium alloy (titanium, aluminum, vanadium), stainless steel and DLC (Diamond Like Carbon)
- In contact with central circulatory system and/or central nervous system

### POSSIBLE COMPLICATIONS

Lack of adequate quantity and/or quality of bone, lack of implant initial stability, excessive bone compression, surgical trauma, infection, generalized diseases and surgical and prosthetic malpractice are some potential causes for failure of Osseointegration. Too early load of an implant, overload from misfitting prosthetics, poor occlusion and articulation, parafunction and trauma may cause loss of an implant. Further complications can include bleeding, chronic pain and peri-implantitis that can lead to implant failure or postoperative infection and inflammation. It is not possible to ensure a 100% Osseointegration success rate.

### ⚠ WARNINGS

- Immediately and early loaded implants can lead to alveolar process deficiencies, recession, contour and peri implant complications
- Narrow diameter implants, platforms, components and angled abutments, are not recommended for single-crown use in the posterior region
- Component and instrument availability, including dimensional aspects, height and angulations, laboratorial and prosthetic possibilities, anatomical, mesiodistal and interocclusal limitations, and specific clinical case requirements, should be verified during planning, before surgical and prosthetic procedures
- Mechanical stress, including originating from abutment and components modification or preparation, may affect the implant and components stability and accuracy
- Take care when selecting the implant length and drills to avoid damage or perforation of structures, e.g. membranes and nerves. Pre-surgical image evaluation, intra oral exams and patient mouth opening verification should be executed during planning, before surgical and prosthetic procedures
- Expansion techniques, grafting, early and immediate loading, smoking, diabetes, anemia, osteoporosis, nutritional deficiencies, long-term treatment with non-steroid anti-inflammatory, hyperparathyroidism, antimetabolic medication, irradiated bone, severe rheumatoid polyarthritis, infection, local or systemic diseases, drug and alcohol abuse, psychoses, compromised immune system, poor health conditions and factors interfering with bone healing increase the chances of implant failure
- Patients identified as at-risk for Creutzfeldt-Jakob disease (CJD) and related infections should be treated with single-use instruments. Therefore, devices that have been in use or suspected of use in patients with CJD after surgery must be disposed of according to current national recommendations.
- The patient must follow dietary, oral hygiene, pre and post medication instructions

## ADVERSE EVENT

If an adverse event occurs (death, significant deterioration of the patient's health status or threat to public health) in relation to the device, the qualified professional shall immediately inform the manufacturer and the local competent authority. To inform the manufacturer, access [www.pibranemark.com](http://www.pibranemark.com) contact area, or download and submit the completed Complaint Registration Form to [support@pibranemark.com](mailto:support@pibranemark.com)

## INFORMATION FOR PATIENT

Before and after surgery, inform patient about contraindications, possible complications, regular follow-up requirements and dietary, oral hygiene, pre-medication and post-medication instructions. Patients should avoid hot beverages immediately following the implant procedure. Smoking, parafunctional habits, local and systemic factors or diseases can increase the chances of implant failure. Should complications occur, such as bleeding, chronic pain, mobility or postoperative infection and inflammation, the patient should consult a qualified professional. Instruct the patient to immediately report any adverse event to the responsible professional, manufacturer and the local competent authority.

## STERILIZATION

P-I Dental Components are for single use and must not reprocessed nor re-sterilized by user. P-I Dental Components identified as sterile are packed sterile in accordance with EN ISO 11137 and have been gamma sterilized. The expiry date must be observed. P-I Dental Components identified as non-sterile must be sterilized before use according to the User Instructions for Processing Non-sterile Devices below. Re-sterilization of single use devices by the user is not permitted and the manufacturer disclaims any responsibility for this, regardless of the sterilization method. Re-use of single use devices is not permitted and can result in cross contamination or impairment of the device performance.

## STORAGE AND TRANSPORTATION

Ensure that the devices are stored and transported under dry conditions, away from sunlight, and that the original packaging is intact.

## GENERAL PRECAUTIONS

The tissue integration implantation method as described as Osseointegration should only be utilized by qualified and experienced professionals trained in this method. Workshops are available for clinicians to learn more about the method. Contact your local distributor for more information. Each patient must be carefully examined and evaluated by the responsible professional to determine psychological, physical and health status to decide on the suitability and indication, or not, of the treatment with dental implants for the specific patient. The P-I Dental Implant Systems have specific design characteristics for combining Implants, Abutments, Components and Instruments. Combining Components and Instruments that are not original P-I devices, not configured or dimensioned for correct combination, modified, not well-maintained or not sterilized when required, can lead to mechanical failure of Components and Implants, damage to tissue, loss of Osseointegration and compromised results. Close cooperation between the implant surgeon, restorative dentist and dental laboratory technician is essential for success.

The P-I Dental Implant Systems may be used and restored prosthetically only by dentists, doctors, oral and maxillofacial surgeons and dental technicians who are familiar with dental implantology, including diagnosis, preoperative planning, surgical technique and prosthetic restoration.

The devices are small and care must be taken to avoid the devices from being swallowed or aspirated by the patient. To secure Instruments, fixate surgical wire or floss in the designated hole where possible.

Instruments with sharp or pointed edges must be handled with care to avoid injury and contamination.

The target population for the medical devices are individuals that have fully completed their growth phase. All contraindications must be observed.

## PROCEDURAL PRECAUTIONS

## PROSTHETIC

To optimize clinical application and minimize adverse reactions and risks for the patient, planning and appropriate implant surgical and prosthetic techniques should be followed.

P-I Instruments should be used for the installation of P-I Dental Components. When any signs of damage, wear, distortion and corrosion, the Instruments must be disposed of. Instruments must be cleaned and sterilized before every use.

Appropriate selection, handling and installation of Components should be adopted to avoid excessive compression and overload to bone, soft tissue, Implant interface and Components. It is important to achieve proper and even stress distribution between the Implants to minimize transverse and excessive loading. Occlusion adjustments, for adequate contact after Osseointegration is achieved, or absence of contact during the healing and provisional phases, should be observed. The opposite jaw may require adjustments. The prosthetic framework must have a passive fit. Components must fit precisely to the Implant or Component interface, without interference of tissue or residues. The verification by

radiographic image can be performed to confirm the correct fitting. If dentures are used after implant installation, they should be alleviated and relined with a soft liner to reduce premature loading.

The removal of Components installed on Implants and Analogs with the MT interface is performed after the removal of the fixation Screw by using the MT Retriever.

Parallel Healing Abutments with height 4.5mm on the MT interface have height marks (1.5, 3 and 4,5mm) to assist in the dimensional selection of subsequent Components.

Angled Conical Abutments are only available for the Conical Abutment 4.1 which has a Platform Ø4.8mm.

The Cylinders Conical Abutment in Titanium can also be used to support prostheses manufactured by a laboratory using CAD/CAM systems. Starting with the first channel, counting from the seating area, the height of the Cylinders Conical Abutment can be adjusted by cutting the channel in the lower section, resulting in Cylinder Conical Abutment in Titanium heights of: 4.5; 5.5; 7, 8 and 9mm. The option of engaging Cylinders Conical Abutment (ENG) is only available for straight Conical Abutments in the AMP and MT interfaces.

In case of blasting, bonding and cementation procedures, residues must be completely removed to avoid contamination of tissues and interference with Components, Screws and Instruments. The use of appropriate dental materials observing the respective manufacturer instructions should be observed. Seating areas of Components and Screws and areas exposed to tissue should not be blasted. Before installing the Components, ensure that screw access and seating areas remain unchanged without modification.

Components should be cleaned and sterilized before use in patients according to cleaning instructions below.

Osseointegration time varies significantly from patient to patient and depends on the specific clinical case aspects and the individual bone healing pattern. The decision as to when to provisionalize and load an Implant, is the sole responsibility of the qualified professional.

The use of a Torque Wrench is recommended to reduce the risk of screws becoming loose. When Components are immediately used on Implants, the torque applied to the Component and Screws must not exceed the initial torque value obtained in the Implant installation. Do not overtighten. Separately supplied Screws should be used in laboratorial steps.

Recommended torques for Components by interface are as follows:

Ncm	⬡ HEX	⊙ AMP	⦿ MT	Driver
Abutments Cylinders over Implant Links	35*	25	25	Ø1.2**
Cylinders - Conical Abutments	15	15	15	Ø1.2
Cover Screws Healings Abutments Healing Caps Impression Copings Scan Bodies	Manual	Manual	Manual	Ø1.2***

\* Except HEX Ø3.5 Components and Angled Conical Abutments = 25 Ncm.

\*\* Except Straight Conical Abutment, Driver Ø2.0.

\*\*\* Except clicked and cemented Abutment Cemented Cylinder Components.

**Melting point:** CoCrMo - Cobalt Chromium Molybdenum | 1,390°C (2,535°F). POM | 165°C (330°F).

Casting and overcasting laboratorial procedures can cause dimensional changes to Components resulting in misfitting of Components and Screws leading to mechanical failure of devices and prosthetics. Verification of Components modified by a laboratory should be thoroughly verified and are the responsibility of the qualified professional.

## USER INSTRUCTIONS FOR PROCESSING NON-STERILE DEVICES

### CLEANING AND DISINFECTION PROCEDURES

All P-I devices identified in the table above as non-sterile must be cleaned and sterilized before use. After cleaning, sterilization by autoclave and drying is recommended. The user should process the Instruments and Components according to EN ISO 17665-1.

## PRECAUTIONS

Cleaning procedures must be performed before clinical application and if debris is visible.

- Do not use warm water
- The devices must be placed in cold tap water (room temperature) for 60 minutes
- When using automated cleaning and disinfection, avoid direct contact of devices to each other
- The devices may not be cleaned using hydrogen peroxide or high chlorine content or containing oxalic acid. Disinfection solution should be aldehyde free
- Do not apply unreasonable force, especially levering and over-bending
- After cleaning and sterilization, Components should be handled only with proper sterile instruments to avoid any contamination

## MANUAL CLEANING PROCEDURES

### MANUAL PRE-CLEANING PROCEDURES

- Rinse the devices under cold tap water until all visible soiling is removed. Firmly adhering soiling should be removed with a soft brush
- Place devices in an alkaline cleaner (e.g. alkaline cleaner 0.5% Neodisher MediClean) for 10 minutes and maximum temperature of 40°C (104°F)
- Rinse the devices under cold tap water to remove the detergent
- Manual drying with a lint-free cloth

### MANUAL DISINFECTION PROCEDURES

- Full immersion of the device in a disinfectant (e.g. Cidex OPA) at 20±2°C (68±3.6°F) for 12 minutes
- Submerge for 1 minute in cold demineralized water
- Extensive flushing with cold demineralized water to remove remaining disinfectants

### AUTOMATIC CLEANING PROCEDURES

- Pre-Cleaning for 4 minutes with cold tap water
- Cleaning with an alkaline cleaner (e.g. alkaline cleaner 0.5% Neodisher MediClean) for 6 minutes and maximum temperature of 55°C (131°F)
- Neutralization with warm deionized water > 40°C (>104°F) for 3 minutes
- Rinse with warm deionized water > 40°C (>104°F) for 2 minutes

### STEAM STERILIZATION

Before sterilization, the original packaging should be removed and the devices should be single-wrapped in sterilization paper. The following sterilization procedures are recommended:

Steam Sterilization Procedure	Parameters
Fractionated Pre-vacuum Method	132°C to 135°C (270°F to 275°F) for at least four minutes. Dry for 20 to 40 minutes.
Gravitation Method	121°C to 123°C (250°F to 254°F) for at least 90 minutes. Dry for 15 to 30 minutes.

Variables that may affect drying times include: loading density of the case/tray, device configuration, total contents of the sterilizer, steam quality, equipment maintenance and others.

Procedures established in national standards and instructions of manufacturers of approved products and equipment for cleaning and sterilization should also be taken in consideration.

### WARNINGS FOR CLEANING PROCEDURES

- Do not sterilize damaged devices
- Check devices for damage or distortion after cleaning and sterilization

## ⚠ PRECAUTIONS

- Avoid contact between devices during the process of cleaning, sterilizing and drying
- Instruments and Components should never be stored without complete drying, in order to avoid oxidation

## 📦 PACKAGING AND HANDLING

All P-I devices identified in the table above as sterile are sterilized by irradiation and marked "sterile". Re-sterilization is prohibited. Components may not be used in the event the outside packaging and/or inside sterile barrier packaging are damaged or signs of previous opening are evident, or if the expiration date has passed.

The sterile barrier should be intact and opened by pulling the tab with the blister facing down.

## 🗑 DISPOSAL

After use, the P-I devices are potentially hazardous to health as they may be contaminated through contact with blood, body fluids, bone or other tissue. The devices and packaging must be treated and disposed of in accordance with local health, environment and governmental legislation and regulations.

Products that have sharp or pointed edges must be handled and disposed of in such a way that there is no danger to third parties.


## 📄 DOCUMENTATION

Each P-I device is traceable using the lot number. The label contains all the corresponding details to transfer to the patient's records.

## 🧲 MAGNETIC RESONANCE IMAGING SAFETY INFORMATION

The P-I Dental Implant Systems have not been evaluated for safety and compatibility in an MRI (Magnetic Resonance Imaging) environment. It has not been tested for heating, migration, or image artifact in the MRI environment. The safety of P-I Dental Implant Systems in the MRI environment is unknown. Scanning a patient who has this device may result in patient injury.

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**SIC invent Sweden AB**