Developed By P-I Brånemark

INSTRUCTIONS FOR USE ENGLISH

P-I Dental Implants

This Instructions for Use applies to all P-I Dental Implants included in the P-I Dental Implant Systems.

Implants	Material	HEX	AMP	MT	Single Use	Sterile
		HEX-S	AMP-S	MT-F		
M+N	Titanium	G4	G4	G4	◦	◦
OSPOL	Titanium	G4	G4	G4	◦	◦

Material

G4 | commercially pure titanium grade 4 (ASTM F67).

The P-I Dental Implants have the following dimensional range: diameter from 3.3 to 5.0mm, platform from 3.3 to 5.1mm, length from 6 to 18mm. Narrow Implants are defined to have platform diameters of 3.3 and 3.5mm. Short Implants are 6 and 7mm in length.

Consult the device label and options at www.pibranemark.com

DISCLAIMER OF LIABILITY

Some conditions, whether combined or not, represent contraindications, limitations and risks (relative and absolute) for the treatment of patients with implants. The procedures for placement of implants are complex and require specialized training which lies in the responsibility of the respective qualified professional. The professional is the sole responsible to determine if the treatment, procedures and devices are suitable for the specific clinical case and patient. P-I disclaims all and any liability and responsibility with regards to the professional decisions and procedures when using the P-I devices.

Some devices may not be available or cleared for clinical use in your region.

INDICATIONS FOR USE

The P-I Dental Implants are endosseous implants part of the P-I Dental Implant Systems and are intended to be used for anchorage of crowns, bridges or dentures to the surrounding bone in the upper and lower jaws to reestablish esthetic and masticatory function. The procedures may range from prosthetic replacement of a single tooth to entire arch of bridgework. Crowns and bridges can be screw-retained and/or cemented to implants and/or abutments.

CONTRAINDICATIONS

Pre-operative patient evaluation is necessary to determine any factors which may put the patient at risk from the implant placement procedure itself, or factors that may affect healing capacities of either the bone or associated soft tissue.

Implants should not be used:

- In patients with a health history that contraindicates surgical intervention
- Where bone is of insufficient quality and/or quantity or not available to produce adequate stability and support
- When there is risk for overload, unfavorable jaw relations and/or parafunction
- Where there is no space to place sufficient number of implants in an optimum position to support biomechanical loads
- Where bone quantity and/or quality is not able to provide adequate implant initial stability
- Where there is a proven allergy or hypersensitivity to titanium, one of the metals in the titanium alloy (titanium, aluminum, vanadium), stainless steel and DLC (Diamond Like Carbon)
- In contact with central circulatory system and/or central nervous system

POSSIBLE COMPLICATIONS

Lack of adequate quantity and/or quality of bone, lack of implant initial stability, excessive bone compression, surgical trauma, infection, generalized diseases and surgical and prosthetic malpractice are some potential causes for failure of Osseointegration. Too early load of an implant, overload from misfitting prosthetics, poor occlusion and articulation, parafunction and trauma may cause loss of an implant. Further complications can include bleeding, chronic pain and peri-implantitis that can lead to implant failure or postoperative infection and inflammation. It is not possible to ensure a 100% Osseointegration success rate.

⚠ WARNINGS

- Immediately and early loaded implants can lead to alveolar process deficiencies, recession, contour and peri implant complications
- Narrow diameter implants, platforms, components and angled abutments, are not recommended for single-crown use in the posterior region
- Short implants should be used in conjunction with regular length implants
- Component and instrument availability, including dimensional aspects, height and angulations, laboratorial and prosthetic possibilities, anatomical, mesiodistal and interocclusal limitations, and specific clinical case requirements, should be verified during planning, before surgical and prosthetic procedures
- Mechanical stress, including originating from abutment and components modification or preparation, may affect the implant and components stability and accuracy
- Take care when selecting the implant length and drills to avoid damage or perforation of structures, e.g. membranes and nerves. Pre-surgical image evaluation, intra oral exams and patient mouth opening verification should be executed during planning, before surgical and prosthetic procedures
- Expansion techniques, grafting, early and immediate loading, smoking, diabetes, anemia, osteoporosis, nutritional deficiencies, long-term treatment with non-steroid anti-inflammatory, hyperparathyroidism, antimitotic medication, irradiated bone, severe rheumatoid polyarthritis, infection, local or systemic diseases, drug and alcohol abuse, psychoses, compromised immune system, poor health conditions and factors interfering with bone healing increase the chances of implant failure
- Patients identified as at-risk for Creutzfeldt-Jakob disease (CJD) and related infections should be treated with single-use instruments. Therefore, devices that have been in use or suspected of use in patients with CJD after surgery must be disposed of according to current national recommendations.
- The patient must follow dietary, oral hygiene, pre and post medication instructions

ADVERSE EVENT

If an adverse event occurs (death, significant deterioration of the patient's health status or threat to public health) in relation to the device, the qualified professional shall immediately inform the manufacturer and the local competent authority. To inform the manufacturer, access www.pibranemark.com contact area, or download and submit the completed Complaint Registration Form to support@pibranemark.com

INFORMATION FOR PATIENT

Before and after surgery, inform patient about contraindications, possible complications, regular follow-up requirements and dietary, oral hygiene, pre-medication and post-medication instructions. Patients should avoid hot beverages immediately following the implant procedure. Smoking, parafunctional habits, local and systemic factors or diseases can increase the chances of implant failure. Should complications occur, such as bleeding, chronic pain, mobility or postoperative infection and inflammation, the patient should consult a qualified professional. Instruct the patient to immediately report any adverse event to the responsible professional, manufacturer and the local competent authority.

STERILIZATION

P-I Dental Implants are for single use and must not reprocessed nor re-sterilized by user. P-I Dental Implants are packed sterile in accordance with EN ISO 11137 and have been gamma sterilized. The expiry date must be observed. Re-sterilization of single use devices by the user is not permitted and the manufacturer disclaims any responsibility for this, regardless of the sterilization method. Re-use of single use devices is not permitted and can result in cross contamination or impairment of the device performance.

⚠ STORAGE AND TRANSPORTATION

Ensure that the devices are stored and transported under dry conditions, away from sunlight, and that the original packaging is intact.

⚠ GENERAL PRECAUTIONS

The tissue integration implantation method as described as Osseointegration should only be utilized by qualified and experienced professionals trained in this method. Workshops are available for clinicians to learn more about the method. Contact your local distributor for more information. Each patient must be carefully examined and evaluated by the responsible professional to determine psychological, physical and health status to decide on the suitability and indication, or not, of the treatment with dental implants for the specific patient. The P-I Dental Implant Systems have specific design characteristics

for combining Implants, Abutments, Components and Instruments. Combining Components and Instruments that are not original P-I devices, not configured or dimensioned for correct combination, modified, not well-maintained or not sterilized when required, can lead to mechanical failure of Components and Implants, damage to tissue, loss of Osseointegration and compromised results. Close cooperation between the implant surgeon, restorative dentist and dental laboratory technician is essential for success.

The P-I Dental Implant Systems may be used and restored prosthetically only by dentists, doctors, oral and maxillofacial surgeons and dental technicians who are familiar with dental implantology, including diagnosis, preoperative planning, surgical techniques and prosthetic restoration.

The devices are small and care must be taken to avoid the devices from being swallowed or aspired by the patient. To secure Instruments, fixate surgical wire or floss in the designated hole where possible.

Instruments with sharp or pointed edges must be handled with care to avoid injury and contamination.

The target population for the medical devices are individuals that have fully completed their growth phase. All contraindications must be observed.

PROCEDURAL PRECAUTIONS

SURGICAL AND PROSTHETIC

To optimize clinical application and minimize adverse reactions and risks for the patient, planning and appropriate implant surgical and prosthetic techniques should be followed.

Recommended Surgical Sequence for healed sites according to Implant geometry and bone tissue density:

HYBRID IMPLANTS: HEX, HEX-S, AMP, AMP-S AND MT

mm	Drill Ø	2.2	2.8	3.4	3.8	4.6	4.8
	rpm	600 - 1,200					
Implant Ø							
3.3		•					
3.75		•					
4.0		•					
4.8		•	•				
5.0		•	•				

= Soft = Normal = Dense

The use of Dense Drills (15 – 50 rpm) should be considered during all surgical preparation and before the Hybrid Implant insertion, regardless of Implant type and bone density, with the objective of not exceeding 50 Ncm of insertion torque value and to avoid excessive cortical bone compression.

INTERPOLATED IMPLANTS: MT-F

mm	Drill Ø	2.2	2.8	3.4	3.8	4.6	
	rpm	600 - 1,200					
Implant Ø							
3.3		•					
3.75		•					
4.1		•					
4.8		•	•				

6 mm Tip = Soft = Normal = Dense

The subsequent Conical Drill, in terms of diameter, should be considered with a drilling depth of 6 mm, in order to not exceed 70 Ncm of insertion torque value. The use of Dense Drills (15 – 50 rpm) can also be considered to lower the insertion torque value.

PROCEDURAL PRECAUTIONS (CONTINUED)

P-I Instruments should be used for surgical preparation and installation of P-I Dental Implants. When any signs of damage, wear, distortion and corrosion, the Instruments must be disposed of. Instruments must be cleaned and sterilized before every use.

It is recommended to apply delicate and coordinated in-and-out axial movements when using Initial and Conical Drills, for conventional and guided surgery, without lateral movements, with irrigation directed to the entrance margin of the Instrument in the surgical site, at low rotation (600 – 1,200 rpm). Osteotomy should be performed at the total depth of the Implant planned position. Drill marks are <1.0mm longer than the nominal measurement. Drills marks in mm:

The removal of cortical bone with Dense Drills should be performed at very low rotation (15 - 50 rpm), without irrigation. Dense Drills can also be considered to gradually prepare surgical sites, bone regularization, widening of the cortical region and in post extraction sites, to avoid cortical bone compression or deviation of Implant position during insertion.

Depending on the clinical situation and surgical procedure, the responsible professional may decide to use the Soft Tissue Punch, in conventional or guided surgery, and Spade Drill and/or Round Burr, in conventional surgery, preceding the steps described in the Surgical Sequence.

Always verify the effective total drilling depth of all Instruments, including the depth resulting from the association of Drills using Stops, before performing osteotomy drilling steps.

In cases of post extraction and sinus lift Implant installation, surgical preparation procedures should be carefully evaluated and determined by the responsible qualified professional.

Care should be taken to prevent soft tissue from being captured by the Implant during installation.

Drills should always be inspected before use to ensure no mechanical damage and should be replaced after approximately 20 uses or when cutting ability diminishes.

Guide Pins and Depth Probes assist in dimensional verification of sites during surgery.

The Implant Insertion Driver is used to take the Implant from the internal packaging to install the Implant in the prepared site, at low rotation, with handpiece and/or Torque Wrench, without irrigation. The Implant Insertion Driver must be completely adapted to the Implant internal hexagonal indexation and, when used with Torque Wrench, the upper hexagonal portion of the Implant Insertion Driver should be entirely connected to the hexagon of the Driver Adapter. The Implant Insertion Driver dots, and its upper hexagon, are indexed to the Implants hexagonal internal indexation. The removal of the Implant Insertion Driver should be carefully performed, axially, without force and lateral movements, to avoid loss of stability of the installed Implant. The Drill Extension may not be used with Implant Insertion Drivers.

Drills and Instruments with handpiece connection are only to be used when handpiece interface is according to ISO 1797.

Do not apply force to connect, remove or correct the Instrument direction during use in conventional and guided surgery.

Appropriate selection, handling and installation of Components should be adopted to avoid excessive compression and overload to bone, soft tissue, Implant interface and Components. It is important to achieve proper and even stress distribution between the Implants to minimize transverse and excessive loading. Minimum distance between two implants or between Implant and teeth should be observed.

Osseointegration time varies significantly from patient to patient and depends on the specific clinical case aspects and the individual bone healing pattern. The decision as to when to provisionalize and load an Implant, is the sole responsibility of the qualified professional.

For more information, see P-I Rotating Instruments and P-I Components Instructions for Use.

Do not exceed the maximum recommended insertion torque value for each Implant stated below:

	HEX	HEX-S	AMP	AMP-S	MT	MT-F
Ncm	50	50	50	50	50	70

PACKAGING AND HANDLING

All P-I devices identified in the table above as sterile are sterilized by irradiation and marked "sterile". Re-sterilization is prohibited. Implants may not be used in the event the outside packaging and/or inside sterile barrier packaging are damaged or signs of previous opening are evident, or if the expiration date has passed.

The sterile barrier should be intact and opened in an aseptic environment by pulling the tab. Inside the sterile barrier packaging, there is the inner Implant container which is opened by rotating the lid counterclockwise. Following the opening of the Implant container, the Implant is taken with the Implant Insertion Driver to be installed in the prepared surgical site. The Implant should be installed directly after opening the inside sterile barrier packaging. Exposure and contact of the Implant with surfaces other than the inner titanium Implant container must be avoided.

DISPOSAL

After use, the P-I devices are potentially hazardous to health as they may be contaminated through contact with blood, body fluids, bone or other tissue. The devices and packaging must be treated and disposed of in accordance with local health, environment and governmental legislation and regulations.

Products that have sharp or pointed edges must be handled and disposed of in such a way that there is no danger to third parties.


DOCUMENTATION

Each P-I device is traceable using the lot number. The label contains all the corresponding details to transfer to the patient's records.

MAGNETIC RESONANCE IMAGING SAFETY INFORMATION

The P-I Dental Implant Systems have not been evaluated for safety and compatibility in an MRI (Magnetic Resonance Imaging) environment. It has not been tested for heating, migration, or image artifact in the MRI environment. The safety of P-I Dental Implant Systems in the MRI environment is unknown. Scanning a patient who has this device may result in patient injury.

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SIC invent Sweden AB