Developed By P-I Brånemark

INSTRUCTIONS FOR USE ENGLISH

P-I Dental Rotating Instruments

This Instructions for Use applies to the below listed P-I Dental Rotating Instruments, non-sterile, included in the P-I Dental Implant Systems.

	Materials		♻️ Reprocessable
Instruments			
Driver	Stainless Steel	Non-sterile	°
Drill	Stainless Steel	Non-sterile	°
Punch	Stainless Steel	Non-sterile	°
Stop	Titanium Grade 5	Non-sterile	°
Extension	Stainless Steel	Non-sterile	°

Some devices may feature Diamond Like Carbon surface.

The P-I Dental Rotating Instruments have dimensional options. Guided Surgery Instruments are identified as "GS".

Consult the device label and options at www.pibranemark.com

DISCLAIMER OF LIABILITY

Some conditions, whether combined or not, represent contraindications, limitations and risks (relative and absolute) for the treatment of patients with implants. The procedures for placement of implants are complex and require specialized training which lies in the responsibility of the respective qualified professional. The professional is the sole responsible to determine if the treatment, procedures and devices are suitable for the specific clinical case and patient. P-I disclaims all and any liability and responsibility with regards to the professional decisions and procedures when using the P-I devices.

Some devices may not be available or cleared for clinical use in your region.

INDICATIONS FOR USE

The P-I Dental Rotating Instruments are part of the P-I Dental Implant Systems and are intended to be used for anchorage of crowns, bridges or dentures to the surrounding bone in in the upper and lower jaws to reestablish esthetic and masticatory function. The procedures may range from prosthetic replacement of a single tooth to entire arch of bridgework. Crowns and bridges can be screw retained and/or cemented to implants and/or abutments.

INDICATIONS

Spade Drill is used for marking the cortical bone and assists on the preliminary definition of osteotomy direction.

Round Burr is used for marking, regularizing and contouring cortical bone.

Soft Tissue Punch is used to remove soft tissue providing access to cortical bone in flapless surgical procedures and can be used through the respective Sleeve positioned in the surgical guide.

Crestal Drill - GS is used for the perforation of crestal bone to facilitate the subsequent osteotomy drilling steps and can be used through the respective Sleeve positioned in the surgical guide.

Initial Drill is used for the initial perforation of crestal bone and can be used to determine the osteotomy depth according to planned implant position. The Initial Drill - GS can be used with Stops and Stops – GS.

Conical Drill is used for the finalization of osteotomy according to bone density and implant dimensions at the planned position. The Conical Drill - GS can be used with Stops and Stops – GS.

Dense Drill is used to remove cortical bone and for the gradual preparation of the osteotomy in the crestal region. It can be used for bone regularization, adjustments and widening of post extraction sites. Dense Drills – GS can be used through the respective Sleeve positioned in the surgical guide.

Pin Drill - GS is used to provide access by drilling the bone through the respective Pin Sleeve for the fixation of surgical guide when using Guided Surgery Pin

Stop is used to limit osteotomy drilling depth in guided surgery. When used through the respective Sleeve installed in the surgical guide, it also provides for a drilling direction.

Drill Extension is used to extend the Drill shaft to provide access to bone when there is limited mesiodistal space to perform the osteotomy drilling steps.

Implant Insertion Driver is used to take the Implant from the internal packaging to install the Implant in the prepared surgical site, at low rotation, with handpiece and/or Torque Wrench. The Implant Insertion Driver – GS is used through the respective Sleeve positioned in the surgical guide.

CONTRAINDICATIONS

Pre-operative patient evaluation is necessary to determine any factors which may put the patient at risk from the implant placement procedure itself, or factors that may affect healing capacities of either the bone or associated soft tissue.

Rotating Instruments, Implants and Components should not be used:

- In patients with a health history that contraindicates surgical intervention
- Where bone is of insufficient quality and/or quantity or not available to produce adequate stability and support
- When there is risk for overload, unfavorable jaw relations and/or parafunction
- Where there is no space to place sufficient number of implants in an optimum position to support biomechanical loads
- Where bone quantity and/or quality is not able to provide adequate implant initial stability
- Where there is a proven allergy or hypersensitivity to titanium, one of the metals in the titanium alloy (titanium, aluminum, vanadium), stainless steel and DLC (Diamond Like Carbon)
- In contact with central circulatory system and/or central nervous system

POSSIBLE COMPLICATIONS

Lack of adequate quantity and/or quality of bone, lack of implant initial stability, excessive bone compression, surgical trauma, infection, generalized diseases and surgical and prosthetic malpractice are some potential causes for failure of Osseointegration. Too early load of an implant, overload from misfitting prosthetics, poor occlusion and articulation, parafunction and trauma may cause loss of an implant. Further complications can include bleeding, chronic pain and peri-implantitis that can lead to implant failure or postoperative infection and inflammation. It is not possible to ensure a 100% Osseointegration success rate.

⚠️ WARNINGS

- Immediately and early loaded implants can lead to alveolar process deficiency, recession, contour and peri implant complications
- Component and instrument availability, including dimensional aspects, height and angulations, laboratorial and prosthetic possibilities, anatomical, mesiodistal and interocclusal limitations, and specific clinical case requirements, should be verified during planning, before surgical and prosthetic procedures
- Take care when selecting the implant length and drills to avoid damage or perforation of structures, e.g. membranes and nerves. Pre-surgical image evaluation, intra oral exams and patient mouth opening verification should be executed during planning before surgical and prosthetic procedures
- Expansion techniques, grafting, early and immediate loading, smoking, diabetes, anemia, osteoporosis, nutritional deficiency, long-term treatment with non-steroid anti-inflammatory hyperparathyroidism, antimetabolic medication, irradiated bone, severe rheumatoid polyarthritis, infection, local or systemic diseases, drug and alcohol abuse, psychoses, compromised immune system, poor health conditions and factors interfering with bone healing increase the chances of implant failure
- Patients identified as at-risk for Creutzfeldt-Jakob disease (CJD) and related infections should be treated with single-use instruments. Therefore, devices that have been in use or suspected of use in patients with CJD after surgery must be disposed of according to current national recommendations.

- The patient must follow dietary, oral hygiene, pre and post medication instructions

ADVERSE EVENT

If an adverse event occurs (death, significant deterioration of the patient's health status or threat to public health) in relation to the device, the qualified professional shall immediately inform the manufacturer and the local competent authority. To inform the manufacturer, access www.pibranemark.com contact area, or download and submit the completed Complaint Registration Form to support@pibranemark.com

INFORMATION FOR PATIENT

Before and after surgery, inform patient about contraindications, possible complications, regular follow-up requirements and dietary, oral hygiene, pre-medication and post-medication instructions.

Patients should avoid hot beverages immediately following the implant procedure. Smoking, parafunctional habits, local and systemic factors or diseases can increase the chances of implant failure. Should complications occur, such as bleeding, chronic pain, mobility or postoperative infection and inflammation the patient should consult a qualified professional. Instruct the patient to immediately report any adverse event to the responsible professional, manufacturer and the local competent authority.

STERILIZATION

P-I Dental Instruments are delivered non-sterile and must be cleaned and sterilized before every use according to the User Instructions for Processing Non-sterile Devices below. Cleaning and sterilization are the responsibility of qualified professional and the manufacturer disclaims any responsibility when not observed which can result in cross contamination or impairment of the device performance.

⚠️ STORAGE AND TRANSPORTATION

Ensure that the devices are stored and transported under dry conditions, away from sunlight, and that the original packaging is intact.

⚠️ GENERAL PRECAUTIONS

The tissue integration implantation method as described as Osseointegration should only be utilized by qualified and experienced professionals trained in this method. Workshops are available for clinicians to learn more about the method. Contact your local distributor for more information. Each patient must be carefully examined and evaluated by the responsible professional to determine psychological, physical and health status to decide on the suitability and indication, or not, of the treatment with dental implants for the specific patient. The P-I Dental Implant Systems have specific design characteristics for combining Implants, Abutments, Components and Instruments. Combining Components and Instruments that are not original P-I devices, not configured or dimensioned for correct combination, modified not well-maintained or not sterilized when required, can lead to mechanical failure of Components and Implants, damage to tissue, loss of Osseointegration and compromised results. Close cooperation between the implant surgeon, restorative dentist and dental laboratory technician is essential for success.

The P-I Dental Implant Systems may be used and restored prosthetically only by dentists, doctors, oral and maxillofacial surgeons and dental technicians who are familiar with dental implantology, including diagnosis, preoperative planning, surgical techniques and prosthetic restoration.

The devices are small and care must be taken to avoid the devices from being swallowed or aspired by the patient. To secure Instruments, fixate surgical wire or floss in the designated hole where possible.

Instruments with sharp or pointed edges must be handled with care to avoid injury and contamination.

The target population for the medical devices are individuals that have fully completed their growth phase. All contraindications must be observed.

PROCEDURAL PRECAUTIONS

SURGICAL AND PROSTHETIC

To optimize clinical application and minimize adverse reactions and risks for the patient, planning and appropriate implant surgical and prosthetic techniques should be followed.

Recommended Surgical Sequence for healed sites according to Implant geometry and bone tissue density:

HYBRID IMPLANTS: HEX, HEX-S, AMP, AMP-S AND MT

mm	Drill Ø	2.2	2.8	3.4	3.8	4.6	4.8
Implant Ø							
3.3			• (S-W)				
3.75			• (S)	(N-D)			
4.0			• (S)		(N-D)		
4.8			•	•	(S)	(N-D)	
5.0			•	•	(S)		(N-D)

rpm 600 - 1,200

S = Soft
 N = Normal
 D = Dense

The use of Dense Drills (15 – 50 rpm) should be considered during all surgical preparation and before

the Hybrid Implant insertion, regardless of Implant type and bone density, with the objective of not exceeding 50 Ncm of insertion torque value and to avoid excessive cortical bone compression.


INTERPOLATED IMPLANTS: MT-F

mm	Drill Ø	2.2	2.8	3.4	3.8	4.6
Implant Ø	rpm	600 - 1,200				
3.3		S ↓	D ↓			
3.75		S		N ↓	D ↓	
4.1		S	N	D		
4.8		S			S	N ↓
	↓ 6 mm Tip	S = Soft	N = Normal	D = Dense		

The subsequent Conical Drill, in terms of diameter, should be considered with a drilling depth of 6 mm, in order to not exceed 70 Ncm of insertion torque value. The use of Dense Drills (15 – 50 rpm) can also be considered to lower the insertion torque value.

P-I Instruments should be used for surgical preparation and installation of P-I Dental Implants. When any signs of damage, wear, distortion and corrosion, the Instruments must be disposed of. Instruments must be cleaned and sterilized before every use.

It is recommended to apply delicate and coordinated in-and-out axial movements when using Initial and Conical Drills, for conventional and guided surgery, without lateral movements, with irrigation directed to the entrance margin of the Instrument in the surgical site, at low rotation (600 – 1,200 rpm). Osteotomy should be performed at the total depth of the Implant planned position. Drill marks are <1.0mm longer than the nominal measurement. Drills marks in mm:

The removal of cortical bone with Dense Drills should be performed at very low rotation (15 - 50 rpm), without irrigation. Dense Drills can also be considered to gradually prepare surgical sites, bone regularization, widening of the cortical region and in post extraction sites, to avoid cortical bone compression or deviation of Implant position during insertion.

Depending on the clinical situation and surgical procedure, the responsible professional may decide to use the Soft Tissue Punch, in conventional or guided surgery, and Spade Drill and/or Round Burr, in conventional surgery, preceding the steps described in the Surgical Sequence.

Always verify the effective total drilling depth of all Instruments, including the depth resulting from the association of Initial and Conical Drills – GS and Stops, in advance of initiating osteotomy drilling steps.

The accuracy, stability, effective offset and dimensional limitations of the planned and laboratory manufactured surgical guide, with or without Sleeves, should be verified before surgery and is the responsibility of the qualified professional.

In cases of post extraction and sinus lift Implant installation, surgical preparation procedures should be carefully evaluated and determined by the responsible qualified professional.

Care should be taken to prevent soft tissue from being captured by the Implant during installation.

Drills should always be inspected before use to ensure no mechanical damage and should be replaced after approximately 20 uses or when cutting ability diminishes.

Guide Pins and Depth Probes assist in dimensional verification of sites during surgery

The Implant Insertion Driver is used to take the Implant from the internal packaging to install the Implant in the prepared site, at low rotation, with handpiece and/or Torque Wrench, without irrigation. The Implant Insertion Driver must be completely adapted to the Implant internal hexagonal indexation and, when used with Torque Wrench, the upper hexagonal portion of the Implant Insertion Driver should be entirely connected to the hexagon of the Driver Adapter. The Implant Insertion Driver dots, and its upper hexagon, are indexed to the Implants hexagonal internal indexation. The removal of the Implant Insertion Driver should be carefully performed, axially, without force and lateral movements, to avoid loss of stability of the installed Implant. The Drill Extension may not be used with Implant Insertion Drivers.

The Implant Insertion Drivers – GS, Narrow and Regular are used in the respective Sleeves – GS (Narrow or Regular) positioned in the surgical guide. The upper lines of the widest diameter portion indicate 9mm and 10.5mm (+1.5mm) surgical guide offset, respectively.

Drills and Instruments with handpiece connection are only to be used when handpiece interface is according to ISO 1797.

Do not apply force to connect, remove or correct the Instrument direction during use in conventional and guided surgery.

Appropriate selection, handling and installation of Components should be adopted to avoid excessive compression and overload to bone, soft tissue, Implant interface and Components. It is important to achieve proper and even stress distribution between the Implants to minimize transverse and excessive loading. Minimum distance between two implants or between Implant and teeth should be observed.

Osseointegration time varies significantly from patient to patient and depends on the specific clinical case aspects and the individual bone healing pattern. The decision as to when to provisionalize and load an Implant, is the sole responsibility of the qualified professional.

For more information see P-I Dental Implants and P-I Dental Components Instructions for Use.

Do not exceed the maximum recommended insertion torque value for each Implant stated below:

	○	⊙	⊙	⊙	⊙	⊙
	HEX	HEX-S	AMP	AMP-S	MT	MT-F
Ncm	50	50	50	50	50	70

Conventional Surgery

Initial and Conical Drills for conventional surgery have a total length of 32mm, and the drilling depth is visually determined by observing the marks as pictured above.

Initial and Conical Drills – GS can be used in conventional surgery by observing the marks as pictured above.

Guided Surgery

Initial and Conical Drills – GS have lengths: Short = 35.5mm and Long = 40mm.

Stops – GS allow for 9mm surgical guide offset and should be used with the respective Initial and Conical Drills – GS (Narrow or Regular, Short or Long).

Stops – GS have a Narrow or Regular (Short or Long) configuration and when associated to the respective Initial and Conical Drill – GS (Narrow or Regular – Short or Long) will orient the direction and limit the drilling depth of the osteotomy through the respective Sleeve – GS (Narrow or Regular) positioned in the surgical guide.

Crestal and Dense Drill - GS (Narrow or Regular) allow for a 9mm surgical guide offset and should be used without applying lateral force or out-of-center movements through the respective Sleeves – GS (Narrow or Regular) positioned in the surgical guide. Crestal Drill – GS is used before the steps described in the Surgical Sequence above.

Pin Drill - GS should be used through the Pin Sleeve – GS positioned in the surgical guide, without applying lateral force. The drilling depth is determined during planning when associating the Pin Sleeve – GS with the the Pin Drill – GS which has a total drilling depth of 13mm when the embedded stop area touches the upper portion of the Pin Sleeve - GS. It is the responsibility of the qualified professional to determine the adequate number of Guided Surgery Pins to stabilize a surgical guide and avoid undesired distortion and movement.

Stops

The Stops should be fully connected to the respective Initial and Conical Drills – GS using manual clockwise rotation until the Stop rests on the provided stop platform. To remove the Stop, apply manual counterclockwise rotation or use the Stop Removal Tool – GS when necessary. The installation and removal of Stops should only be performed when the Drills are not connected to the handpiece.

USER INSTRUCTIONS FOR REPROCESSABLE DEVICES

CLEANING AND DISINFECTION PROCEDURES

All P-I devices identified in the table above as reprocessable must be cleaned and sterilized before use. After cleaning, sterilization by autoclave and drying is recommended. The user should reprocess the Instruments according to EN ISO 17665-1.

Immediately after use, the Instruments must be fully disassembled, placed in cold water (room temperature) or disinfection solution, avoid any drying on of blood, tissue or any other secretion. For manual cleaning, rinse the device under cold tap water until all visible soiling is removed. Firmly adhering soiling should be removed with a soft brush. After disinfection and drying, reassemble the Instruments and pack them each in a suitable sterilization paper.

PRECAUTIONS

- By using machined cleaning and disinfection, avoid direct contact of the instruments to each other
- Cleaning procedures must begin within 1 hour after the surgical application
- Immediately after use, place the Instruments and Components in cold water or disinfection solution. Do not use warm water or detergent as this may lead to fixation of the residues on the device, which can affect the outcome of the subsequent cleaning step
- For stainless steel instruments, always use solvents appropriate for this material
- The devices may not be cleaned using hydrogen peroxide or high chlorine content or containing oxalic acid. Disinfection solution should be aldehyde free
- Do not apply unreasonable force, especially leveraging and over-bending
- Do not mix steel instruments and stainless steel instruments on any cleaning, disinfection and sterilization process

MANUAL PRE-CLEANING PROCEDURES

- The devices must be placed in cold tap water (room temperature) for 60 minutes

MANUAL CLEANING PROCEDURES

- Rinse the devices under cold tap water until all visible soiling is removed. Firmly adhering soiling should be removed with a soft brush
- Place devices in an alkaline cleaner (e.g. alkaline cleaner 0.5% Neodisher MediClean) for 10 minutes and maximum temperature of 40°C (104°F)
- Rinse the devices under cold tap water to remove the detergent
- Manual drying with a lint-free cloth

MANUAL DISINFECTION PROCEDURES

- Full immersion of the device in a disinfectant (e.g. Cidex OPA) at 20±2°C (68±3.6°F) for 12 minutes
- Submerge for 1 minute in cold demineralized water
- Extensive flushing with cold demineralized water to remove remaining disinfectants

AUTOMATIC CLEANING PROCEDURES

- Pre-Cleaning for 4 minutes with cold tap water
- Cleaning with an alkaline cleaner (e.g. alkaline cleaner 0.5% Neodisher MediClean) for 6 minutes and maximum temperature of 55°C (131°F)
- Neutralization with warm deionized water > 40°C (>104°F) for 3 minutes
- Rinse with warm deionized water > 40°C (>104°F) for 2 minutes

STEAM STERILIZATION

Instruments and Components delivered in a non-sterile condition must be cleaned, disinfected, and sterilized before the initial and each subsequent use. Before sterilization, the original packaging should be removed and the devices should be single-wrapped in sterilization paper. The following sterilization procedures are recommended:

Steam Sterilization Procedure	Parameters
Fractionated Pre-vacuum Method	132°C to 135°C (270°F to 275°F) for at least four minutes. Dry for 20 to 40 minutes.
Gravitation Method	121°C to 123°C (250°F to 254°F) for at least 90 minutes. Dry for 15 to 30 minutes.

Variables that may affect drying times include: loading density of the case/tray, instrument configuration, total contents of the sterilizer, steam quality, equipment maintenance and others.

Procedures established in national standards and instructions of manufacturers of approved products and equipment for cleaning and sterilization should also be taken in consideration.

⚠ WARNINGS FOR CLEANING PROCEDURES

- Do not sterilize corroded or rusty instruments
- Check instruments for corrosion after sterilization

PRECAUTIONS

- Avoid contact between instruments during the process of cleaning, sterilizing and drying
- Do not apply temperatures on stainless steel instruments above 135°C (275°F) on any operation
- Do not mix steel instruments and stainless steel instruments on any cleaning, disinfection and sterilization process
- Instruments and Components should never be stored without complete drying, in order to avoid oxidation

⚠ PACKAGING AND HANDLING

All P-I Instruments identified in the table above are non-sterile and must be cleaned and sterilized before each use. Instruments may not be used in the event the packaging is damaged or signs of previous opening are evident.

The packaging should be intact and opened by pulling the tab with the blister facing down.

DISPOSAL

After use, the P-I devices are potentially hazardous to health as they may be contaminated through contact with blood, body fluid, bone or other tissue. The devices and packaging must be treated and disposed of in accordance with local health, environment and governmental legislation and regulations.

Products that have sharp or pointed edges must be handled and disposed of in such a way that there is no danger to third parties.


DOCUMENTATION

Each P-I device is traceable using the lot number. The label contains all the corresponding details to transfer to the patient's records.

MAGNETIC RESONANCE IMAGING SAFETY INFORMATION

The P-I Dental Implant Systems have not been evaluated for safety and compatibility in an MRI (Magnetic Resonance Imaging) environment. It has not been tested for heating, migration, or image artifact in the MRI environment. The safety of P-I Dental Implant Systems in the MRI environment is unknown. Scanning a patient who has this device may result in patient injury.

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SIC invent Sweden AB