

Kit Composition



Developed By P-I Brånemark

	Description	Advanced	Start-up	Specialist Kit	Specialist Kit	Prosthetic
Application	MT	▲	Ø3.75	Ø3.75	Ø4.0	▲
	AMP-S and AMP	▲	Ø3.75 4.0	Ø3.75	Ø4.0	▲
	HEX-S and HEX	▲	Ø3.75 4.0	Ø3.75	Ø4.0	▲

Code		Description	Advanced	Start-up	Specialist Kit	Specialist Kit	Prosthetic
131139	Drivers	Implant Insertion Driver - Medium	•	•	•	•	
131140		Implant Insertion Driver - Long	•	•			
131141		Implant Insertion Driver - HEX 3.5 Medium	•				
131142		Implant Insertion Driver - HEX 3.5 - Long	•				
131146		Implant Insertion Driver - Medium - Manual			•	•	
141138	Drills	First Drill - 2.2	•	•	•	•	
141146		Conical Drill - 2.8	•	•	•	•	
141148		Conical Drill - 3.4	•	•	•		
141314		Conical Drill - 3.8	•	•		•	
141152		Conical Drill - 4.6	•				
141315		Conical Drill - 4.8	•				
141213		Dense Drill - 3.3	•				
141316		Dense Drill - 3.75	•	•	•		
141215		Dense Drill - 4.0	•	•		•	
141317		Dense Drill - 4.8 5.0	•				
131010	Instruments	Hexagonal Driver - Short - 1.2					•
131011		Hexagonal Driver - Medium - 1.2	•	•			•
131012		Hexagonal Driver - Long - 1.2					•
131016		Conical Abutment Driver - Short					•
131017		Conical Abutment Driver - Medium	•	•			•
131131		Retriever - MT	•	•			•
131136	Accessories	Torque Wrench	•	•			
131028		Drill Extension	•	•			
131114		Guide Pin 2.2 2.8	••	••			
131115		Guide Pin 2.8 3.8	••	••			
131117		Stainless Steel Tray	a	a			
131134		Polymer Tray	a	a			
131138		Compact Tray	a	a	•	•	•
Total items (including 1 tray)			24	18	7	7	8

KIT COMPOSITION-007-ENG-2018-07-19

a - Select Preference